Raypa®

Avgda. del Vallès, 322
Pol. Ind. "Els Bellots"
08227 Terrassa
(Barcelona) España

Tel. +34.93.783.07.20
Tel. +34.93.783.03.44
Fax +34.93.731.37.17

raypa@raypa.com
www.raypa.com

2015

Raypa®

R. Espinar, S.L.

CATÁLOGO GENERAL

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Raypa

● **LABORATORY EQUIPMENT
AND INSTRUMENTATION**

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R. Espinar, S.L. is a spanish company dedicated to manufacture Instrumentation and Equipment for Laboratory since 1970. All this time has become our best argument to guarantee quality and durability at a reasonable price. Our Sales network covers more than 60 countries all over the world, supporting our aim of exporting company.

Our facilities in Terrassa are equipped with the suitable technology to develop, manufacture and complete assembling of the products included in our catalogue.

All of our products have **CE** mark. Our manufacture process accomplishes with regulations 2006/95/CE (Electric Safety), 2004/108/CE (Electromagnetic Compatibility), and 97/23/CE (Pressure Containers), externally certified by TÜV Rheinland number 1027.

In addition to the above mentioned, and to guarantee a global quality service, **R. Espinar, S.L.** is an ISO 9001 certified company.

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
Due to the continuous improvement process, information herein is subject to change without notice. **R. Espinar, s.l.** assumes no liability for any errors that may appear in this document.

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Research & Development • Cutting-edge technology

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NORMATIVAS / Regulations / Directives

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Our instruments fulfill the following requirements of the european directives:

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Regulations

- EN 61010-1 Safety electrical equipment requirements.
- EN 61010-2-010 Particular requirements for heating equipment.
- EN 61010-2-041 Particular requirements for autoclaves using steam.
- EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring.
- EN 61326 Electromagnetic compatibility.

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According to the test results, the manufacturer is entitled to mark the pressure equipment produced within the range of this module with the mark:

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*have been audited and approved by Quality Management Systems
in compliance with the requirements of the standard:*

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
Communication software RS232 for data recording.

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Most of Raypa instrumentation is prepared to be connected with a computer. Software basic functions are:

- Data collection: data are captured by the computer to carry out reports, worksheet processing, etc.

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CERTIFICADOS / Certificates / Certifications

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Calibration certificates. ENAC standards traceability. Protocols for OQ, IQ validations.

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Intended use

* AES/AHS-N

Suitable to sterilize: culture media, glass, liquids, plastic, metallic elements, waste bags, etc.

** AE-DRY/AHS-DRY

Suitable to sterilize: unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, liquids, metallic objects, waste bags.

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Intended use

N Class

- Suitable to sterilize: culture media, glass, liquids, plastic, metallic elements, waste bags, etc.

S Class

- Suitable to sterilize: unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, metallic objects, waste bags, etc. AH-21 L is suitable to sterilize liquids.

B Class

- Suitable to sterilize: unwrapped and wrapped solids, porous objects, and hollow ware.

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
MODELOS AES / Models AES / Modèles AES

- The sterilization with saturated water steam is the most common used method. The Raypa® autoclaves use three different methods to purge and remove the air from the chamber:

1) **Purge by gravity displacement.** The air is moved away when water saturated steam is produced.

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
MODELOS AE-DRY / Models AE-DRY / Modèles AE-DRY

2) **Purge by vacuum.** The air is pumped out using a vacuum pump and replaced by saturated steam.

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
CLASE B / B Class / Classe B

3) **Purge by fractionated vacuum.** The air is pumped out with several vacuum pulses to make sure that air has been completely removed and replaced by water steam injected from a steam generator (B class).

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Advantages

- Small and easy use.
- Programmable auto-start (99 h).
- Connection to PC.
- Control software (optional).

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Intended use


- **Suitable to sterilize:** culture media, glass, liquids, plastic, metallic elements, waste bags, etc.

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Features

- External case and top made of AISI-304 stainless steel.
- Chamber and lid made of 18/10 stainless steel.
- Temperature regulation by microprocessor.
- Includes RS232 connection for PC.
- Manual valves to drain and to unsteam.

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Safety

- Safety valve and thermostat.
- Open lid blocking system.

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Supplied with the following accessories

- Heater cover tray and stainless steel basket.

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CV-8
CESTO INCLUIDO / Basket included / Panier compris

REFERENCIA Reference Référence	CAPACIDAD l Capacity l Capacité l	Ø × H cm	H × A × F cm H × W × D cm H × L × P cm	TEMP. °C	PRESIÓN MÁX. bar Max. pressure bar Pression max. bar	W	kg	RANGO TIEMPO min. Time range min. Rang de temps min.
AES-8	8	22 × 18	43 × 34 × 34	100-127	1,5	1000	11	0-500

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
Features

- External case, top, chamber and lid made of AISI-304 stainless steel.
- Total control by microprocessor: ten sterilization programs. Six programs are user free.
- Includes RS232 port to connect a printer or to a PC.
- It includes inlet in chamber for an external validation probe.
- Supplied with protecting grid for the heaters. Unsteaming pushbutton. Manual valve to drain.

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Safety

- Safety valve and thermostat.
- Hydraulic door blocking system while existing positive pressure in the chamber.
- Open lid sensor.

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Models AES-12, AES-28, AES-50 and AES-75 are supplied with 230 V mono-phase 230 V (mono-phase 110V optional). Models AES-110 and AES-150 are supplied with 400 V three-phase supply (three-phase 230 V optional).

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Intended use

- **Suitable to sterilize:** unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, liquids, metallic objects, waste bags.

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
Features

- External case, top, chamber and lid made of AISI-304 stainless steel.
- Total control by microprocessor: ten sterilization programs. Six programs are user free.
- Feeding water reservoir. Automatic water filling.
- Includes RS232 port to connect a printer or to a PC.
- It includes inlet in chamber for an external validation probe.
- Supplied with protecting grid for the heaters.

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Safety

- Safety valve and thermostat.
- Hydraulic door blocking system while existing positive pressure in the chamber.
- Safety pressure switch.
- Open lid sensor.

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Models AE-28 DRY, AE-50 DRY and AE-75 DRY are supplied with 230 V mono-phase 230 V (mono-phase 110V optional). Models AE-110 DRY and AE-150 DRY are supplied with 400 V three-phase supply (three-phase 230 V optional).

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Advantages

- Automatic water feeding.
- Steam generator.
- Fractionated vacuum.
- Connection to PC.
- Control software (optional).
- Built in printer (optional).

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Intended use

- **B Class - Suitable to sterilize:** unwrapped and wrapped solids, porous objects, and hollow ware.

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
Features

- Total control by microprocessor. Test cycles: Vacuum Test, Bowie & Dick Test / Helix Test.
- Clean and used water reservoir.
- Ten predefined programs and four user free.
- Three PT100 control probes.
- Includes RS232 port to connect to a printer/PC.

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Safety

- Safety valve, thermostat and safety pressure switch.
- Hydraulic door blocking system.
- Open lid sensor.

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- Tray support for five trays ④.

- AH-21-B: Three stainless steel perforated trays ref. BAH-21 ③.
- AHS-50/75-B: Two stainless steel perforated trays ref. BAH-50/75.
- A holding clamp ⑤.
- Flexible tube with quick connection to drain condensation reservoir ②.
- Auxiliary plastic tray ①.

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
Intended use

- **S Class - Suitable to sterilize:** unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, metallic objects, waste bags. AH-21 L is suitable to sterilize liquids.

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Features

- Pre-vacuum prior sterilization and vacuum system for final drying.
- Controlled by microprocessor.
- 4.5 l and 2 litres feeding and condensates reservoirs. Automatic chamber filling.
- It includes RS232 connection for printer or PC.

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Safety

- Safety valve and thermostat.
- Door hydraulic blocking system.
- Safety pressure switch.
- Open lid sensor.

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
- Tray support for four trays ④.
- Three stainless steel perforated trays: ref. BAH-21 ③.
- A holding clamp ⑤.
- Flexible tube with quick connection to drain condensation reservoir ⑥.
- Bag support ①.
- Auxiliary tray ②.

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
Intended use

- **S Class - Suitable to sterilize:** unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, metallic objects, waste bags, etc.

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Features

- Pre-vacuum prior sterilization and vacuum system for final drying.
- Controlled by microprocessor.
- Feeding and condensates reservoirs. Capacity: 4.5 litres. Chamber automatic filling.
- It includes RS232 connection for external printer.

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Safety

- Safety valve and thermostat.
- Door hydraulic blocking system.
- Open lid sensor.

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
Advantages


- Selection for solids and liquids.
- Connection to PC.
- Control software (optional).
- Programmable start up (500 h).

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Intended use

- **N Class - Suitable to sterilize:** culture media, glass, liquids, plastic, metallic elements, waste bags, etc.

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Features

- Controlled by microprocessor.
- 7.5 litres feeding water reservoir.
- It includes RS232 connection to PC.

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Safety

- Safety valve and thermostat, door blocking system while positive pressure exists in the chamber.
- Open lid sensor.

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- Tray support for four trays ②.
- Three stainless steel perforated trays: ref. BAH-21 ③.
- A holding clamp ④.
- Flexible tube with quick connection to drain condensation reservoir ⑤.
- Auxiliary tray ①.

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Features

- Chamber, lid and closure made of AISI-304 stainless steel.
- Total control by microprocessor: ten sterilization programs. Six programs are user free.
- Includes RS232 port to connect to a printer/PC.

- Unsteam push button.
- Manual drain valve.
- 7 litres feeding water reservoir.
- Manual filling of the vessel.

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
- Support for five trays ②.
- Two perforated trays made of stainless steel ref. BAH-75 / BAH-50 ③.
- Auxiliary tray ①.

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Intended use


- Suitable to sterilize: unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, liquids, metallic objects, waste bags.

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Features

- Chamber, lid and closure made of AISI-304 stainless steel.
- Total control by microprocessor: ten sterilization programs. Six programs are user free.
- Automatic water filling.
- Automatic drying.
- Includes RS232 port to connect to a printer/PC.

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/ Stainless steel Baskets. Not perforated / *Paniers en acier inoxydable. Non perforés*

REFERENCIA Reference Référence	PARA AUTOCLAVE For autoclave Pour autoclave	DIMENSIONES mm Dimensions mm Dimensions mm	N.º CESTOS AUTOCLAVE No. baskets autoclave Nº paniers autoclave	PESO kg Weight kg Poids kg
CCI-12	AES-12	Ø 225 × 215	1	1
CCI-28	AES-28/AE-28 DRY	Ø 280 × 200	2	1,5
CCI-28	AES-50/AE-50 DRY	Ø 280 × 200	3	1,5
CCI-75	AES-75/AE-75 DRY	Ø 380 × 260	2	2,5
CCI-75	AES-110/AE-110 DRY	Ø 380 × 260	3	2,5
CCI-150	AES-150/AE-150 DRY	Ø 475 × 360	2	4

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Flexible "heart" probe PT-100

- Recommended for sterilizing liquids in flasks of 1 litre or bigger.
- The control is done from the microprocessor.
- This accessory should be installed in factory.

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Referencia: IT

Thermal paper Printer

- It prints program number, cycle number, temperature, date and hour of the run.
- Error messages are also printed.
- To be installed in our facilities.

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Referencia: SOFT-RS232

PC connection kit for RS232

- Kit for the connection and use of the communication software between the equipment and the PC, through the RS232 interface.
- It includes: Female USB-RS232 converter. CD with corresponding drivers included.
- Male-male DB-9 cable to connect converter and the equipment.

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
Referencia: ITS

Table top printer

- It prints program number, cycle number, temperature, date and hour of the run.
- Error messages are also printed.

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Referencia: TEST-BD

Bowie & Dick test

- Class B indicator.
- Check the absence of residual air.
- Box of 30 tests.

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Referencia: TEST-HT

Test de Helix

- To check air removal (steam penetration) capability of type B steam steriliser's.

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Kjeldahl

- The Kjeldahl method is a chemical analysis process to determine the contents of nitrogen in a specific sample. It is widely used as a reference method to estimate the protein contents in food.

Soxhlet

- The Soxhlet method is used to determine fat contents in food and some other products. The extraction is carried out through a solvent extraction.

Weende

- The Weende method is commonly used to quantify crude fibre in a food sample after an acidic and alkaline hydrolysis treatment.

Van Soest

- The Van Soest method is used to determine acid detergent raw fibre (ADF) and neutral detergent raw fibre (NDF) contents in a sample.

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/ Equipment for food and water analysis

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Advantages

- Programmable auto-start.
- Microprocessor built-in the case. Nine user free programs. All programs can have up to 18 steps for temperature/ incubation.

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
Ver página 43

Intended use

- Analysis of proteic nitrogen in organic substances as: nourishing products, fertilizers, seeds, beer, meat, sauces, water, petroleum, etc.

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Features

- Quick and safe manipulation of samples.
- Used together with "Scrubber" does not need fume hood.
- Integrated support for sample tubes enabling an easier cooling.
- Sample size solids: 5 g (macro), 1 g (micro). Liquids: 15 ml (macro), 3 ml (micro).

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- One tube holder rack with digestion tubes (tubes included) ①.
- One fume collector ②.

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• Advantages

- Vacuum pump included.
- Saving of refrigerating water in digestion heating blocks.
- Fumes neutralization closed system.
- Filtration and neutralization of fumes.

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Features

- The gas produced during the digestion is extracted by the action of a vacuum pump. First, the gas undergoes a phase of condensation that acts as preliminary extractor for steams and dragged liquids, avoiding a warming or increase of the volume in later washing solution.
- The acid or alkaline steam is washed and neutralized in the next step. In this neutralization phase, most of the unwanted particles are retained by means of granules of activated coal. The used air is taken directly to an extractor or outdoors.
- The neutralization unit "Scrubber", is connected by a flexible tube to the digestion fume collector.

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Advantages

- 20 user programs.
- Pumps calibration by end-user.
- Electric steam generator with water level control.
- Door open sensor.

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Features

- Control of the system by microprocessor with LCD screen.
- Twenty user free programs.
- Additionally, the system has two fixed programs:
 - Ammonium sulphate test.
 - Preheating of the system.

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Safety

- Closed door detection.
- Digestion tube detection.
- Over-temperature thermostat.

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Performance

- Measurement range: de 0.1 to 200 mg of N.
- Nitrogen recovery: $\geq 99.5\%$.
- Reproducibility: $\pm 1\%$.
- Distillation time: 4 min for 100 ml.
- Cooling water consume: 1.2 l/min.
- Electric power: 2000 W (100%).
- Number of programs: 20 user programs.

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Programmable parameters and values

- Water for dilution: 0-200 ml.
- NaOH solution: 0-200 ml.
- Reaction time (delay): 0-30 min.
- Distillation time: 0-30 min or indefinite.
- Steam power regulation: 30-100%.

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- 1 tube series macro Ø 42 x 300 mm.
- 2 reservoirs of 10 litres for H₂O and NaOH.
- Several connection tubes.

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Advantages

- Optional kit for automatic titration.
- Sample suction.
- 20 user programs.
- Pumps calibration by end-user.
- Electric steam generator with water level control.
- Boric acid pump.
- Door open sensor.

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Features

- Control of the system by microprocessor with LCD screen.
- Twenty user free programs.
- Additionally, the system has three fixed programs
 - Ammonium sulphate test.
 - Preheating of the system.
 - Washing.

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Performance

- Measurement range: de 0.1 of 200 mg de N.
- Nitrogen recovery: $\geq 99.5\%$.
- Reproducibility: $\pm 1\%$.
- Distillation time: 4 min for 100 ml.
- Cooling water consume: 1.2 l/min.
- Electric power: 2000 W (100%).
- Number of programs: 20 user programs.

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Programmable parameters and values

- Water for dilution: 0-200 ml.
- NaOH solution: 0-200 ml.
- H_3BO_3 solution: 0-200 ml.
- Reaction time (delay): 0-30 min.
- Distillation time: 0-30 min or indefinite.
- Steam power regulation: 30-100%.

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- 1 tube series macro $\varnothing 42 \times 300$ mm.
- 2 containers of 10 litres for H_2O , NaOH and 1 container of 5 litres for H_3BO_3 .
- Several connection tubes.

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Advantages

- Automated system.
- Results in N or total protein.
- Sample suction.
- 20 user programs.
- Pumps calibration by end-user.
- Electric steam generator with water level control.
- Built in magnetic stirrer for titration.
- Door open sensor.

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Features

- Control of the system by microprocessor with LCD screen.
- Automatic titration of distilled sample.
- Additionally, the system has three fixed programs:
 - Ammonium sulphate test.
 - Preheating of the system.
 - Washing.
- Automatic titrator included.

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Components supplied with the equipment

- Automatic titrator.
- 1 tube series macro $\varnothing 42 \times 300$ mm.
- 2 containers of 10 litres for H_2O , NaOH and 1 container of 5 litres for H_3BO_3 .
- Several connection tubes.

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It includes: magnetic stirrer, reaction vessel and complete titrator with the necessary accessories.

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Titration features

- 5 titration programs.
- pH calibration.
- Burette autocalibration.
- Results in N or total protein.
- Additional applications.

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Titration specifications

- Resolution: 0.01 pH, 1 mV, 0.1 °C.
- Reproducibility: ± 0.01 pH.
- Temperature compensation with probe Cat. Pt100.

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/ Distillation tubes and glass spare parts

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Intended use

- Determination of the volumetric alcoholic degree in wines, musts, beers and spirit drinks of high graduation. Volatile acidity and sorbic acid. According to regulation CEE 2676/90 y CEE 2870/2000.

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Features

- Distillation unit by steam dragging.
- It has a steam generator automatically fed with distilled water, having electronic level control and automatic stop in the case of lack of water in the container.

- Safety thermostat.
- Automatic stop when reaching 200 ml of distillate.
- Manual levelling of the distillate to 200 ml avoiding adding distilled water.
- Automatic stop of refrigerating water at the end of distillation.
- Optical display of minimal level of water in the steam generator.
- Button to add distillate manually.

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Technical data

- Refrigeration water flow: 75 l/h.
- Water consumption (steam generator): 1.6 l/h.
- Distilled water reservoir: 10 l.
- Approx. distillation time: 7 min./200 ml.

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Components supplied with the equipment

- 1 sample tube macro Ø 42 x 300 mm.
- 1 sample tube macro Ø 80 x 300 mm.
- 1 volumetric flask of 200 ml.
- 1 reservoir of 10 litres for H₂O.
- Connection tubing.

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Advantages

- Independent programming for the plate temperature, extraction time "boiling", and extraction time "rinsing".

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Intended use

- Extraction system with solvents, according to Soxhlet and Randall methods for fat and soluble materials determination, in food, feeds, soils, detergents, polymers, paper pulp, fibers, rubbers, textiles, petrochemical products, pharmaceutical products, etc.

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Features

- Furniture in stainless steel, epoxy resin painted.
- Heating by an electric plate with armoured heater to assure homogeneous distribution of heat.
- Control of temperature by microprocessor with P.I.D. action and probe Pt100 class A.
- Over temperature safety control by independent thermocouple.

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Technical Data

- Temperature range: room temp. +5 °C to +250 °C.
- Solvent volume: up to 50 ml.
- Max. sample volume: 25 ml.
- Solvent recovery: 60-70%.
- Reproducibility: ±1%.
- Extraction thimbles: Øext 26 x 60 mm.
- "boiling" extraction time: 0-500 hours.
- "rinsing" extraction time: 0-99 hours.
- Extraction time: 30-60 minutes.
- Electric protection IP65.

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Components supplied with the equipment

- 18 extraction vessels + magnetic adapter.
- Box containing 25 thimbles.
- 7 thimble supports.
- Extraction vessel tong.
- 3 racks for extraction cartridges.
- 3 racks for extraction vessels.
- Aligning rack.
- Insertion handle.

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/ Fibre extraction system Fibertest

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Features

- Extraction and filtration without sample transfer (it avoids lose of sample).
- High reproducibility of conditions and results.
- Up to 36 analysis/day (Weende).
- Versatile use: samples can be dried/weighted at every phase if needed.
- System for six simultaneous samples.
- All electric and mechanic elements are conveniently protected.
- Furniture in stainless steel epoxi painted.

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Technical Data

- Range of measurement: 0.1-100%.
- Reproducibility: $\pm 1\%$ for a fibre contents level of 1-30%.

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Components supplied with the equipment

- Supplied with 6 glass crucibles porosity P-2, test tube crucibles rack, heater lid and handling support.
- Crucible vessel tong.

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Components supplied with the equipment

- Water jet pump, ref. 5700
- 6 glass crucibles, ref. CR-P2.
- Solvent recuperator, ref. RS-EF6.

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Advantages

- Programmable auto-start.
- Microprocessor built-in the case.
- Nine user free programs.
- Each program has 18 different.
- Temperature/incubation time segments.

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Intended use


- C.O.D. (Chemical Oxygen Demand) analysis in waters.
See page 43.

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Features

- Made of aluminium. Housing made for digestion tubes Ø 42 (C.O.D.).
- Heat insulation with ceramic fibre and a double air chamber.
- Circuit protection against corrosive agents.

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Temperature control

- The control has 9 user free programs. Each program has a waiting time to initializing. All programs can have 18 different temperature/incubation time segments. Each segment has a ramp with specific slope and an incubation time.
- Safety system against breakage of temperature sensor.

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Components supplied with the equipment

- 1 test tube rack with support.
- 1 support for cooling tubes.
- 1 set of tubes for C.O.D. with conical inlet 29/32.
- 1 set of cooling tubes for air with conical inlet 29/32.

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

Advantages

- Two programmable speeds.
- Two programmable running times.
- Two handles for transportation.

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Intended use

- Jar test. Water analysis.
- Leaching test.

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Features

- Digital flocculator for water analysis. 6 places.
- Speed control by microprocessor. Adjustable between 15 and 300 rpm.
- Time also programmable between 0 and 99.0 h.

- LCD screen displaying target and real time values.
- Adjustable two speeds and two times.
- Stirring rods adjustable in height.
- Quiet.
- For 1 l vessels (tall or short model).
- Screen with light with two positions for an easier sample checking.

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Advantages

- Car cigarette lighter connection.
- Handy and light.
- With screen light.

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Features

- Portable flocculator for water analysis. 4 places.
- Analogue speed control.

- Adjustable between 20 and 200 rpm.
- Timer 0-60 min.
- Stirring rods adjustable in height.
- Quiet.
- Screen with light for an easier sample checking.
- Car cigarette lighter connection cable.

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Features

- Table top compact water distiller, made of borosilicate glass.
- Mounted on a metallic case, painted with epoxy resin.
- Stainless steel heating element (INCOLOY), protected against over temperature with safety thermostat.
- Production of high purity distilled water, pyrogen free and low conductivity.
- Water production: 4 l/h.
- Safety: norm EN 61010-1.

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/ Reservoir with filling level. Automatic distiller switch off when it's full

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Features

- Table top compact water distiller made of borosilicate glass. It can be installed on the wall.
- Mounted on a metallic case, painted with epoxy resin.

- Quartz heating element with screw for an easy and quick replacement.
- Safety system that disconnects heating in the case of failure in the refrigeration water supply. Automatic re-start when water supply is recovered.
- High purity distilled water production, free of metallic ions and pyrogenic substances.
- Production: 4 l/h.
- Safety: norm EN 61010-1.

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Advantages

- It takes little space.
- Cartridge easily replaced.
- Conductivity gauge in μS .

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Intended use

- Production of deionised water for use in the laboratory.

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Features

- Wall-unit for connection to tap water, suitable for the production of purified and deionised water in small volumes.

- During deionization all inorganic substances are removed, so that an absolute ionic purity is achieved.
- It includes a conductivity meter for actual reading of the produced water conductivity ($\mu S/cm$).

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Advantages

- Without refrigeration water consumption.

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Intended use

- Production of distilled water for use in the laboratory.

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Features

- Air forced distillation.
- Distillation volume: 1.4~ 1.5 l/h.
- Made of thermo-resistant plastic.
- Inner made of stainless steel.
- Capacity: 4 l.
- Over temperature safety.
- Power supply: 220-240 V/110-120 V.
- Power consumption: 750 W.

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Intended use

- Water analysis, B.O.D. test (Biological Oxygen Demand), enzymatic reactions and bacteria culture at low temperature, botanics, aging, cosmetics, samples storage, etc.

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Features

- Adjustable temperatures: from 0 °C to 45 °C.
- Time and temperature controlled by microprocessor.
- Temperature probe: Pt 100 class A.
- Cooling by hermetic compressor.

- Temperature homogeneity: $\pm 1,5$ °C.
- Temperature stability: $\pm 0,15$ °C.
- Two free sockets inside.
- Forced air convection.
- RS232 connection to connect to a PC.

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Intended use

- Preservation, storage and freezing of products.
- Adjustable temperatures to -30°C .
- Temperature homogeneity: $\pm 1.5^{\circ}\text{C}$.

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Features

- Epoxy resin oven painted external case.
- Inner chamber in rough aluminium.
- Inner rails for baskets.
- Folding door.

- Hermetic compressor with anti-vibration system.
- Silent running.
- Control panel with main switch and temperature adjustment by electronic thermostat with digital reading.
- Inner light when opening the lid.
- Cleaning drainage.

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Intended use

- Liquid cooling in thermostatic baths.

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Features

- Table-top equipment made of stainless steel.
- Constant cooling without temperature control.

- Silent hermetic compressor, with air cooled condenser and immersion cooling coil in stainless steel AISI-316.
- Cooling coil dimensions: Ø 32 mm x 160 mm long.
- Refrigerating immersion tube length: 900 mm long.
- Equipment designed to cool down liquids in thermostatic baths, needing a lower temperature than ambient.
- To obtain a real temperature control it must be used together with an immersion thermostat.

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CUBETAS DOBLE PARED / Double wall tanks / Cuves métalliques double paroi

CUBETA DOBLE PARED EN ACERO INOXIDABLE

/ Double wall tank made of stainless steel

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/ Double wall tank made of stainless steel.

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/ INCUBATION AND DRYING
/ INCUBATION ET SÉCHAGE

● INCUBACIÓN Y SECADO / Incubation and drying / Incubation et

● TABLA RESUMEN ESTUFAS INCUBACIÓN

/ Comparative summary culture incubators

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Advantages (serie ID)

- Connection to PC.
- Control software (optional).

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Features

- Natural air convection.
- Adjustable temperature from room temperature +5 °C to +80 °C.
- Temperature regulation by hydraulic thermostat (series I, analogue).
- Temperature regulation by microprocessor (series ID, digital).
- Adjustable safety thermostat.

- Delay & incubation time programmable between 0 and 500 h (series ID, digital).
- Homogeneity and stability according to DIN 12880.
- Safety according to EN 61010-1, EN 61010-2-010.
- Heating by electric elements, assembled in independent chamber, allowing an optimal stabilization of temperature.
- External case epoxy resin painted.
- Inner chamber, inner side door and trays made of stainless steel AISI-304.
- Inner glass door.
- Adjustable aeration device.

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Supplied with

- Two perforated trays made of stainless steel AISI-304 adjustable in height, with anti-tilting system.

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Referencia: SOFT-RS232.

PC connection kit for RS232

- Kit for the connection and use of the communication software between the equipment and the PC.
- It includes: female USB-RS232 converter. CD with corresponding drivers included.
- Male-male DB-9 cable to connect converter with the equipment.
- Installation CD with the software.

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Features

- Natural air convection.
- Adjustable from room temperature +5 °C to +80 °C.
- Temperature regulation by microprocessor.

- Adjustable safety thermostat.
- Homogeneity and stability according to DIN 12880.
- Delay & incubation time programmable between 0 and 500 h.
- Safety according to EN 61010-1, EN 61010-2-010.
- External case epoxy resin painted.
- Inner chamber, inner side door and trays made of stainless steel AISI-304.
- Outer double door.
- Inner double glass door with magnetic closure.
- Adjustable aeration device.

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Supplied with

- Two perforated trays made of stainless steel AISI-304 adjustable in height. With anti-tilting system.

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PC connection kit for RS232

- Kit for the connection and use of the communication software between the equipment and the PC.
 - It includes: female USB-RS232 Converter. CD with corresponding drivers included.
 - Male-male DB-9 cable to connect converter with the equipment.
 - Installation CD with the software.
- Reference: SOFT-RS232.**

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Features

- Natural air convection.
- Temperature adjustable by hydraulic thermostat from 40 °C to 250 °C - DO series, analogue.
- Temperature adjustable by microprocessor from room + 5 °C to 250 °C - DOD series, digital
- Adjustable safety thermostat

- Delay & incubation time programmable between 0 and 500 h (series DOD, digital).
- Homogeneity and stability according to DIN-12880.
- Safety according to EN-61010-1, EN-61010-2-010.
- Heating by electric elements, assembled in independent chamber, allowing an optimal stabilization of temperature.
- External case epoxy resin painted.
- Inner chamber, inner side door and trays made of stainless steel AISI-304.
- Silicone sealing gasket.
- Adjustable aeration device.

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Referencia: SOFT-RS232.

PC connection kit for RS232

- Communication software between the equipment and the PC.
- It includes: female USB-RS232 Converter. CD with corresponding drivers included.
- Male-male DB-9 cable to connect converter with the equipment.
- Installation CD with the software.

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
Features

- Forced air convection.
- Adjustable temperature from room temperature +5 °C to +250 °C.
- Temperature regulation by microprocessor.
- Adjustable safety thermostat.

- Delay & incubation time programmable between 0 and 500 h.
- Temperature homogeneity and stability according to DIN 12880.
- Safety according to EN 61010-1, EN 61010-2-010.
- Heating by electric elements, assembled in independent chamber, allowing an optimal stabilization of temperature.
- External case epoxy resin painted.
- Inner chamber, inner side door and trays made of stainless steel AISI-304.
- Silicone sealing gasket.
- Adjustable aeration device.

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Features

- Adjustable temperatures: from +35 °C to +200 °C.
- Temperature homogeneity: ± 4 °C.
- Temperature stability: ± 1 °C.

- Temperature display accuracy: ± 2 °C.
- Maximum admissible vacuum: 0,01 mbar.
- Safety according to: EN 61010-1, EN 61010-2-010.
- Double safety glass against implosions according to regulations.
- Safety thermostat: class 2 incorporated from factory.
- Time and temperature controller: digital by microprocessor.
- Temperature probe: Pt 100 class A.
- Epoxy resin oven painted external case.
- Interior made in polished stainless steel AISI-304.
- Two perforated extractable shelves made of aluminium with galvanic treatment.

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Advantages

- Overtemp safety thermostat.
- Removable tray. Easy cleaning.

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Intended use

- Incubations, suspension cell culture.

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Features

- Adjustable temperature from room temperature +5 °C to +100 °C.
- Resolution: 0,1 °C.
- Programmable delayed start (0-99,5 h).

- Programmable incubation time (0-99,5 h).
- Safety thermostat with heating system blockage in the case of over-temperature (DIN 12879 class 2).
- Overload protection.
- Inner tank and external case made of stainless steel.
- Shaking width: 40 mm.
- Adjustable shaking speed: 10-150 rpm.
- Adjustable temperature slope °C/min.

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/ Analogue and digital thermostatic baths

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Features

- Series BAE and BOE:
Analogue regulation by hydraulic thermostat with graduated scale in °C and over-temperature warning pilot. Safety thermostat class 2 included.

- Series BAD and BOD:
Temperature regulation by microprocessor with digital read out and P.I.D. action. PT100 probe.
Programmable delayed start (0-500 h).
Programmable incubation time (0-500 h).

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Common features

- Safety under regulation EN 61010-1.
- Safety thermostat with manual rearms according to DIN 12877 class 2.

- Inner tank made of AISI-304 stainless steel.
- External case made of AISI-304 Stainless steel heating elements inside the tank.
- Incorporated drainage.
- Includes heater protecting grid.

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● BAÑOS TERMOSTÁTICOS / Thermostatic baths / Bains thermostatiques

● ACCESORIOS SERIES BAE, BOE, BAD, BOD

/ Accessories Series BAE, BOE, BAD, BOD / Accessoires Séries BAE, BOE, BAD, BOD

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
Features

- Adjustable temperature from room temperature +5 °C to +99.9 °C (or +200 °C).
- Resolution: 0.1 °C.

- Regulation and digital reading of temperature by microprocessor with action P.I.D. and temperature probe PT100 DIN class A.
- Recirculation system.
- Programmable delayed start (0-500 h).
- Programmable incubation time (0-500 h).
- Safety thermostat with heating system blockage in the case of over-temperature (DIN 12879 class 2).
- Inner tank and external case made of stainless steel.
- Includes roof type lid.

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Features

- Inner tank made of AISI-304 stainless steel.
- External case made of AISI-304 stainless steel.
- Frontal control panel.
- Hydraulic thermostat with graduated scale in °C, and pilot light.
- Safety thermostat with manual rearm. It switches off the instrument in the case of over-temperature or low level of liquid in the tank, in accordance with DIN 12877 class 2 regulation.
- Heater protecting grid.

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Features

- Adjustable temperatures on plate surface up to 400 °C.
- Temperature regulation by energy impulses regulator.

- Steel hot plate with surrounding heating elements.
- Stainless steel upper case.
- Control panel with pilot-lights indicator of connection to the mains and heating functioning.
- Removable tray of stainless steel, which permits to use the apparatus as a hot plate.

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Features

- Inner 2l tank and external case made of AISI-304 stainless steel.
- Hydraulic thermostat with graduated scale in °C, and pilot light.

- Safety thermostat with manual rearm. It switches off the instrument in the case of over-temperature or lack of liquid in the tank, in accordance with DIN 12877 class 2 regulation.
- Heater protecting grid.

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● BAÑOS TERMOSTÁTICOS / Thermostatic baths / Bains thermostatiques

● ACCESORIOS SERIES BBA, BBO

/ Accessories series BBA, BBO / Accessoires séries BBA, BBO

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
Advantages

- Economy: low consume of electric energy and detergent.
- Time saving.
- Cavities and spaces with difficult access are cleaned removing all dirty.
- Materials protection: the surface will never be damaged.

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Intended use

- Cleaning of:
Dental and Surgical instruments and materials.
/ Plastic, glass, porcelain, glass materials. / Machinery tools and instruments. / Jewelry, costume jewelry, cutlery, watch industry. / Electronics, printed circuits.
/ Optics, contact lenses. / Automobile accessories, injectors.

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Common features

- Generator of electric energy at high frequency completely transistorized.
- Working frequency 35 kHz.
- Adjustable temperature up to 80 °C on UCI series and up to 90 °C on UCD series.
- Cleaning tray made of stainless steel 18/10.
- External case made of stainless steel AISI-304.
- Heating element with attached elements to the tray.
- Complete-half wave selector. It allows less power consume in some applications.
- Draining valve incorporated.
- Series UCI: analogue regulation of temperature by hydraulic thermostat.
- Series UCD: temperature regulation by microprocessor with digital read out.

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● ACCESORIOS SERIES UCI, UCD

/ Accessories series UCI, UCD / Accessoires séries UCI, UCD

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Features

- Electronic thermostat controlled by microprocessor with digital reading and selection of the temperature. An adjustable safety thermostat blocks the heating system in case of over-temperature (DIN 12879 class 2).
- Temperature sensor PT100 class A.
- Waiting time for a delayed start is programmable from 0 to 500 h.

- Temperature maintenance can be programmed from 0 to 500 hours.
- Re-circulation pump.
- External case made of stainless steel AISI-304.
- All the elements in contact with the liquid are made of stainless steel AISI-304.
- Fixing by a bosshead.
- Pump pressure: 130 mbar.

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Features

- Electronic thermostat controlled by microprocessor with digital reading and selection of the temperature. An adjustable safety thermostat blocks the heating system in case of over-temperature (DIN 12879 class 2).
- Temperature sensor PT100 class A.

- Waiting time for a delayed start is programmable from 0 to 500 h.
- Temperature maintenance can be programmed from 0 to 500 h.
- Circulation pump with external circulation capability.
- Cooling coil to connect to water supply and maintain the ambient temperatures.
- Tank capacity: 8 litres.
- Cover opening 200 x 130 mm.
- Pump pressure: 130 mbar.

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Features

- Electronic thermostat controlled by microprocessor with digital reading and selection of the temperature. An adjustable safety thermostat blocks the heating system in case of over-temperature (DIN 12879 class 2).
- 1/8 HP compressor.

- Temperature sensor PT100 class A.
- Waiting time for a delayed start is programmable from 0 to 500 h.
- Temperature maintenance can be programmed from 0 to 500 h.
- Circulation pump and external circulation.
- Tank capacity: 8 litres.
- Cover opening 200 × 130 mm.
- Pump pressure: 130 mbar.

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Features

- For vessels up to 15 litres. Sturdy, for continuous use.
- Moulded case painted with epoxy.

- Protected against spillage IP54.
- Support for aluminium rod.
- Diameter plate: Ø135 mm.
- Multi point speed control (up to 1500 rpm) and heating (model AG-15 B).
- Connection for electronic regulation (model AG-15 B).

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Advantages

- Waiting time and timer up to 99 h.
- Quick heating (AG-15 V).

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Features

- For vessels up to 15 litres. Sturdy, for continuous use.
- Moulded case painted with epoxy.
- Protected against spillage IP54.
- Support for aluminium rod.

- Diameter plate: 135 mm (AG-15 E).
- Vitro ceramic plate 160 x 160 mm in model AG-15 V (no corrosion).
- Speed control (up to 1500 rpm) and heating control by microprocessor.
- Connection for external probe.
- Adjustable temperature slope °C/min.

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
Advantages

- Proximity sensor.
- Microprocessor control.

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Features

- Epoxy resin oven painted external case IP54.
- Start up and stop switch.

- Sensor/continuous switch to work by proximity sensor or continuously.
- Single tube mixing cup for mixing one tube at a time.
- Speed regulation knob.
- Orbital width: 3 mm.

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
Features

- Silent. For continuous use.
- Speed regulation by microprocessor. LCD display. to 300 r.p.m.
- Progressive start and stop.
- Orbit width: 20 mm.
- Max. load: 15 kg.
- Overload protection.
- Programmable shaking time (0-99 h).
- External metallic case oven painted with epoxy resin.
- Several compatible trays for different applications.

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
Features

- Designed to secure any flat-bottomed glass-ware or plastic-ware: micro or deep well plates, Erlenmeyer or Fernbach flasks, round or square bottles.
- Speeds of up to 300 rpm depending only on the weight of the container.
- Easy to clean.

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Features

- Adjustable inclination of the rotation axis.
- Axis device allows the tube removing while the motors is switched on.
- Rotor not included (see accessories next page).

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

Features

- Adjustable temperatures from 50 °C up to 250 °C.
- For a quick sterilisation of heat resistant materials.
- Natural air convection.

- Epoxy resin oven painted case.
- Tank, cover and perforated tray with handle made of AISI-304 stainless steel.
- Heating by stainless steel electric element.
- Hydraulic thermostat temperature regulation.
- Sterilizing time regulation from 0 to 120 minutes.
- Thermostat and start up pilot-lights.

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Intended use

- Suitable to sterilize needles, rings, sowing grips, etc.
- It avoids possible infections or contaminations by contact and splashing.

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Features

- Quartz incinerator tube, with electric element and thermally isolated with ceramic fibre.
- External grid made of stainless steel for protection.
- Epoxy resin oven painted base. Start up pilot-light.

- Sterilization obtained by means of infrared rays.
- The organic material is incinerated by introducing it into a incandescent quartz tube, avoiding splashing of infectious microorganisms.
- A complete sterilization is reached between 5/8 seconds at 850 °C.

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Intended use

- Furnace to incinerate products, metallurgical assays, etc.

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Features

- Adjustable temperature to 1150 °C.
- Stability of temperature: ± 1 °C.
- Homogeneity of temperature: ± 5 °C.
- The system has a digital control for temperature/time by microprocessor.

- Control P.I.D. with thermocouple type K.
- Nine user free programs.
- Each one has a waiting time to start. 18 different segments can be edited and modified in all user free programs.
- Every segment has a slope and a maintenance time.
- Heating by Ni-Cr thermo-ceramic plates.
- Double body door.
- Fibre-ceramic thermo isolation with low thermo conductivity and low density material.
- Safety system against temperature thermocouple breakage.
- Asbestos free

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MANTAS CALEFACTORAS PARA MATRACES

/ Heating mantles for round flasks / Chauffe-ballons

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Features

- Adjustable temperature up to 450 °C on the surface.
- AISI-304 stainless steel cylindrical external case.
- Insulated glass fibre heater with dimensions according to the corresponding round flask.
- Temperature regulation by energy impulses and pilot-light.

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Features

- Temperature on the surface up to 450 °C. with heating power selector.
- Enamel external case, temperature and corrosion resistant, with heating power selector.
- Sturdy steel hotplate with built in electric heater.

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Features

- Temperature on the surface up to 400 °C with heating power selector.
- Enamel external case, temperature and corrosion resistant, with heating power selector.
- Sturdy steel hotplate with build in electric heater.

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Features

- Adjustable temperature on the plate up to 400 °C.
- Made of AISI-304 stainless steel.

- Hot plates made of rigid steel with homogeneous heating element.
- Independent temperature regulation in every plate by energy impulses controller. Pilot light.
- Adjustable support (height and slope) to hold suitable refrigerating vessels, etc.

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Features

- Temperature on the surface up to 130 °C with heating power selector.
- Enamel external case, temperature and corrosion resistant, with heating power selector.
- Sturdy steel hotplate with build in electric heater.

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Intended use

- Raypa viscometers, models RP1 and RP2, are rotational viscometers for the quick determination of viscosity as specified in ISO 2555 and other ASTM norms.
- Instrument offers a precise digital torque measurement and viscosity measurement results which are 100% compatible to the Brookfield method.

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Features

- Interface RS232 for data transfer to PC or printer, temperature probe PT100, digital display for direct reading of various parameters, software Viscosoft Basic for easy data evaluation, warning of under- or over-range condition, complete spindle set and 12 months warranty.
- Additionally model RP2 offers bi-directional interface RS232 for RP2 control and data evaluation through a PC and the software Viscosoft Plus for the automated control of viscometer and the issuing of rheological studies.

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Compatibility

- ISO 2555, ISO 1652, BS 6075, BS 5350 and ASTM 115-3716 standards define a viscometer in torque, speed and spindle design. Our Raypa viscometers comply with such specifications and are therefore 100% compatible to Brookfield method and able to perform measurements as per such standards.

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Additional information

- Our viscometers are supplied as complete equipment in a suitable carrying case. Standard accessories for each basic equipment: complete set of spindles with storage rack, interconnecting RS232 cable, spindle guard and instructions manual, are included in case.

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Displayed parameters

- Selected speed: rpm.
- Selected spindle: S.P.
- Dynamic viscosity: mPas or cP.
- % Scale (Brookfield): %.
- Sample temperature: °C or °F.
- Auto range to display viscosity limits: mPas or cP.
- Shear ratio S.R. (with special spindles. Only in V2 version): 1/s.
- Shear stress S.S. (with special spindles. Only in V2 version): N/m² or dyne/cm².

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

Referencia: APM.

Features

- Accessory to allow accurate viscosity measurements of small sample volumes (from 8 to 13 ml).
- Works with all models.
- The sample chamber fits into a circulating chamber with controlled temperature. Using a suitable external thermostatic unit, sample can be thermostated between -10°C and $+100^{\circ}\text{C}$.
- It is supplied with a set of 3 or 4 spindles depending on version:
- Version L: TL5 - TL6 - TL7.
- Version R: TR8 - TR9 - TR10 - TR11.

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Referencia: ALV.

Features

- The "Adapter for Low Viscosity Materials" is necessary to obtain accurate measurements on low viscosity materials, both newtonian and no newtonian.
- It can be used with V1 and V2.
- It is supplied with a precision cylindrical spindle (LCP).
- Reproducibility of viscosity measurements from 0.45 mPas.
- Sample volume between 16-18 ml.
- The "Adapter for Low Viscosity Materials" has a circulating chamber containing the sample tube. Through an external element (bath circulator), temperature in sample can be adjusted and controlled (temperature range between -10 °C and +100 °C).

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Referencia: AHM.

Features

- The adapter for helicoidal movement together with special T-shape spindles allows measurement of viscosity in substances as cream, gelatine, paste, wax, mastic, etc. Standard spindles cannot be used.
- The head automatic movement up and down avoids breaking of the sample structure to analyze, obtaining viscosity reliable readings.
- Height adjustable.
- Compatible with all models in our range of viscometers.
- It is supplied with: six suitable spindles with T-shape (PA, PB, PC, PD, PE, PF).

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Features

- Pipettors "Raypipette" have a robust design, suitable for a daily use at the Laboratory.
- They are accurate, variable volume and truly autoclavable (121 °C 20'). There is no need to re-calibrate them after the sterilization process.
- They are equipped with a built-in tip ejector.
- Fully adjustable, selected volume is clearly visible in a digital indicator.
- **Certified according the norms ISO 8655 and DIN 12650.**
- Small volume models are suitable for their use in Molecular Biology micro-techniques.
- Easy maintenance.
- Calibration possible by the end-user.

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● INSTRUMENTAL, SOPORTES Y GRADILLAS

/ Instruments, supports and test tube stands

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/ Bosshead with nickel steel screwing. Maximum opening: Ø 16 mm

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CUBETAS RECTANGULARES CON ESQUINAS REDONDEADAS

/ Rectangular trays with rounded corners

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* Volume divided in 4 with an inner cross.

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Safety base not included (Reference 1-00005).

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SOPORTES PIE PLATO

/ Stand base support / *Support plateau*

REFERENCIA Reference Référence	DIMENSIONES PLATO mm Stand base dim. mm Dimensions plateau mm	DIMENSIONES BARRA mm Rod dimensions mm Dimensions tige mm
2249	220 × 140	8 × 490
2578	250 × 160	12 × 780
27100	270 × 180	15 × 1000

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/ Water jet pump (made of metallic nickel-plated brass)

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Alphabetic index

Alcohol distiller, Enodest	58	Incinerator, bacteria incinerator BAT-850	131
Autoclave, AE	22	Incubators, analogue (models Incuterm)	84
Autoclave, AES	20	Incubators, digital (models Digit ID)	84
Autoclave, AES-8	18	Incubators, digital (models Incuterm-Macro)	86
Autoclave, Clinoclav B-plus	24	Incubators, refrigerated (models IRE)	76
Autoclave, Clinoclav-L	26	Iron gauzes	163
Autoclave, Clinoclav N2	30	Labjack	149
Autoclave, Clinoclav S2-dry	28	Magnetic stirrers	119
Autoclave, Horizontal AHS N	32	Medical graduated cups with handle	153
Autoclave, Horizontal AHS dry	34	Melting spoon	152
Automatic pipettes	147	Micro-spatula spoon/rigid	149
Bars, magnet stirring bars	155	Micro-spatula double rigid	150
Baskets, stainless steel	37	Mohr forceps	160
Baths, analogue (models BAE)	100	Mortar with iron cast and polished	157
Baths, analogue (models BOE)	100	Needle Kolle holder, accessories	155
Baths, thermostatic baths banks (BBA/BBO)	107	Nitrogen distillers, DNP-1500 MP	50
Baths, digital (models BAD)	100	Nitrogen distillers, DNP-2000 MP	52
Baths, digital (models BOD)	100	Nitrogen distillers, DNP-3000	54
Baths, digital with recirculation	103	Ovens, analogue ovens (models DO)	88
Bath, round thermostatic	105	Ovens, digital (models Digit DOD)	88
Bath, shaking BSH	98	Ovens, digital (models DAF)	90
Baths, ultrasonic	109	Ovens, digital (models Digioven-Macro)	92
Beaker tongs	160	Ovens, vacuum oven (model Vacuterm)	94
Blowpipe	161	Petri dishes	152
Bosshead	148	Racks, test tube draining racks	149
Boxes with cover and handle	151	Rectangular trays with round corners	150
Bunsen burner	156	Refrigerated thermostatic unit	114
Burnerhead	156	Retort stand tripod	162
Circulation thermostatic unit	113	Ring with bosshead	151
Clamp, Automatic burette clamp	157	Ring without bosshead	151
Clamp for burettes and tubes	158	Round trays	148
Clamps for burettes with bosshead	158	Round trays with round upper border	148
Cork borer	163	Sand thermostatic bath	106
Crucibles, stainless steel	153	Scissors, metallic	163
Crucible tongs with bow	158	Scrubber	48
Curved forceps	158	Shaker, orbital	124
Curved plastic tongs	159	Sharpener for cork borer	150
Digestion block, MBC-MBCM	44	Sodium spoons	148
Digestion block, (COD) MBDQO	66	Spatula spoon/rigid	153
Digital thermostat DIGI	112	Special supports	161
Dissecting forceps	159	Stainless steel sterilizing boxes	151
Dissecting needles	150	Stand base support	162
Double spatula	150	Suspension hook	154
Fat extractor, Soxtest	60	Teclu burner	156
Fibre extractor, Fibertest	62	Three-prong clamp	160
Flexible spatules with wooden handle	149	Triangle silimanite tubes	164
Flocculator	68	Tripod supports	162
Flocculator, portable	69	Tube rotator	126
Freezers	77	Tube stands, stainless steel wire	154
Funnel with handle	149	Tube stands, stainless steel "z" type	155
Furnaces	132	Tube tongs	161
Fusion spoons	148	Universal clamp	161
Glass tube cutter	153	Vibration spatula	154
Heating mantles	133	Viscometer	140
Hoffmann forceps	160	Vortex mixer	122
Hot air sterilizer	130	Water deionizer DEM-10	73
Hotplates, round PCC	134	Water distiller DEM-4	72
Hotplates, round hotplates banks	136	Water distiller DES-3	71
Hotplates, rectangular hotplates PL	135	Water distiller DES-4	70
Hotplates, rectangular hotplates PR	137	Water jet pumps	164
Immersion refrigerator	78	Weighing spoon	157

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